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EasyLight

Technical parameters - bulb Mitochondriak® 3.0:

Dimensions (mm):	120*120*120	Flicker:	Zero flicker
Power Consumption (Watt):	54W	Flicker Pct:	1,5 – 5 %
LED diode pc:	18	Flicker Index:	0 – 0,03
Weight (kg):	0,38	Material:	Aluminum/Plastic
Angle of enlightenment:	30 °	Lifespan of LED:	50 000 hours
Wavelengths:	630, 670, 760, 810, 830, 850 (1:1:1:1:1:1)	Holder/clamp:	Included in the package
Suitable for use:	Local therapy	Bulb socket:	E27
Input voltage:	AC 110-240 V	Frequency:	50/60 Hz
EMF:	≤ 0,2 mG at the recommended distance	Warranty:	2 years

Bulb intensity:

15 cm: 180 mW/cm²

30 cm: 85 mW/cm²

Warning: you will not see some LEDs lit. Why?

When the bulb is turned on, you will only see about 5 LEDs lit (depending on your eyes perception). The reason why is simple. The human eye can only see light from about 380 to 760 nm. We cannot see UV or infrared. NIR (near infrared light) is not visible to the human eye and therefore you cannot see 760, 810, 830 or 850 nm. However, on our website you will find many measurements with spectrometers where you can see that the LEDs are lit. :)

How to test for photosensitivity?


Before using the any red light therapy panel or bulb, we recommend that you check your photosensitivity. To do this, do the following: Expose an area such as your forearm, 15 cm from the device, to the light for 3 minutes. If you experience any redness or discomfort, you may be photosensitized and we recommend that you consult your doctor or healthcare provider before using the device.

How to install and start a light bulb:

- Remove the protective cover from the bulb.
- Attach the clamp firmly to the table and tighten.
- Attach the gooseneck by screwing it into the socket (from either side).
- Then screw the socket onto the gooseneck.
- Finally, screw the bulb into the socket.
- Once the bulb is installed, you can connect the power cord and turn the switch ON.

How to use a light bulb:

- We recommend wearing safety glasses or not shining the light into your eyes.
- In general, for optimal results, we recommend using the **light at a distance of 30 cm from the body and illuminating for about 10 to 15 minutes per body part. For deeper body parts (bone, organs, etc.) use for 15 to 30 minutes.**
- For the best results, shine the light directly on the skin.
- Consistency is the key to results, so we recommend illuminating 1-2 times a day, at least 5 times a week.

Red light (RED) 630, 670, 760 nm:

Most of the red light is absorbed by your skin, leading to increased collagen production and optimal skin health.

Red light 670nm also best stimulates ATP production and changes the viscosity of water in mitochondria.

Near infrared light (NIR) 810, 830, 850 nm:

Near infrared light penetrates deeper into your body, improving muscle recovery and reducing joint pain. It also stimulates mitochondria and cellular respiration.

General benefits of therapy:

Red light (RED):

- Supports liver function
- Supports circulation
- Beneficial for bones
- Beneficial for kidneys
- Skin renewal
- Supports muscle healing
- Helps blood regeneration
- Improves ATP production

Near infrared light (NIR):

- Accelerates wound healing
- Helps in pain management
- Slows down skin damage
- Strengthens cellular respiration and mitochondrial function
- Relieves scar tissue
- Helps improve metabolism
- Strengthens enzymes
- Improves water viscosity in the body

Warnings and Cautions:

- **EasyLight Mitochondriak® products**, intended for light therapy, **are not medical devices**. The use of our products does not in any way replace a visit to a doctor or healthcare. If you have any health problems, we recommend consulting a specialist.
- Do not look directly into the LEDs or their reflection in the mirror, as this may cause temporary eye irritation.
- Do not use the lamp near water.
- Do not use the lamp if you are photosensitive.
- Do not allow children to use the lamp without supervision.
- We do not recommend using it for pregnant women.
- Do not touch the lamp after therapy and if the lamp overheats, turn it off immediately.
- The lamp must not be used by sleeping or unconscious persons.
- Do not use the lamp on open wounds, as the heat may cause bleeding.
- Do not use the bulb in combination with ointments and balms that contain heat-generating ingredients, as this may cause skin burns.
- Using the bulb in any way other than as described in this manual will void the warranty and may result in serious injury.

Warranty conditions:

We guarantee that this equipment will operate within its specifications for a period of 2 years from the date of delivery.

What is and is not covered by this warranty?

- Two-year warranty
- Damage or wear and tear is not covered. Repairs may be subject to a fee.
- Do not disassemble the device. This will void the warranty.

easylight.sro@gmail.com

www.mitochondriak.com

[instagram/mitochondriak.europe](https://www.instagram.com/mitochondriak.europe)